

PAL Workstation:
Balance-PAL



- For accurate dilutions
- Reagent additions
- Monitoring of weight change reactions (outgassing, digestions)
- Vial transport (various vial sizes)
- Two-way communication Balance-PAL
- Most cost effective; low price tag
- Control it with LEAP Shell

Balance-PAL



- Combines liquid transfers with gravimetric determinations
- Adaptable time schedule
- Adaptable to various vial sizes (2, 5, 10 or 20ml vials)
- Results from weighing sent to Mettler balance link software
- Use automatic barcode scanner
- Small footprint
- High vial capacity
- Works with capped vials and reagent reservoirs
- Add reagent holder and inkwell bottles
- Cooling option to keep samples stable or heating for faster reaction
- Challenge us with your application!



LEAP Technologies PO Box 969 Carrboro, NC 27510 Tel 800 229-8814

www.leaptec.com info@leaptec.com

