

LEAP L-MARK® and CTC X-Type Syringes

(make sure you have appropriate syringe adaptor!)

IF YOU DO NOT SEE AN APPROPRIATE SYRINGE FOR YOUR APPLICATION, CONTACT THE PRODUCT MANAGER

1. GC Applications (GC PAL, COMBI PAL)

Part number	Market	Description	Unit	List Price
LMK.2106302	GC	5ul L-MARK Syringe FN, 26g, AS, fitted plunger	pc	60.00
LMK.2106416	GC	10ul L-MARK Syringe FN, 26g AS fitted plunger	pc	30.00
LMK.2106416-2P	GC	10ul L-MARK Syringe FN, 26g AS fitted plunger, 2/pack	pack	56.00
LMK.2106416-6P	GC	10ul L-MARK Syringe FN, 26g AS fitted plunger, 6/pack	pack	164.00
LMK.2106426	GC	10ul, L-MARK Syringe FN, 23g-26g, dual gauge AS fitted plunger	pc	60.00
LMK.2606403	GC	10ul L-MARK Syringe FN, GT, PTFE, 23g AS	pc	60.00
LMK.2606407	GC	10ul, L-MARK Syringe FN, GT, PTFE, 26g, AS	pc	58.00
LMK.2620706	GC	100ul L-MARK Syringe, FN, GT, PTFE, 26gs, ptAS	pc	93.00
LMK.2620706-2P	GC	100ul L-MARK Syringe, FN, GT, PTFE, 26gs, ptAS 2/pk	pack	177.00
LMK.2620706-6P	GC	100ul L-MARK Syringe, FN, GT, PTFE, 26gs, ptAS 6/pk	pack	505.00
LMK.2620822	GC	250ul L-MARK Syringe, FN, GT, PTFE, 23g, ptAS	pc	119.00
LMK.2620803	GC	250ul L-MARK Syringe, FN, GT, PTFE, 26g, ptAS	pc	127.00
LMK.2620903	GC	500ul L-MARK Syringe, FN, GT, PTFE, 26g, ptAS	pc	127.00
Headspace Syringes				
LMK.2620019	GC	1.0ml L-MARK Syringe FN, GT, PTFE, 22g, pt5, Headspace	pc	120.00
LMK.2620023	GC	2.5ml L-MARK Syringe FN, GT, PTFE, 22g, pt5, Headspace	pc	125.00
LMK.2620032	GC	2.5ml L-MARK Syringe FN, GT, PTFE, 23g, pt5, Headspace	pc	118.00
LMK.2620034	GC	2.5ml L-MARK Syringe FN, GT, PTFE, 26g, pt5, Headspace	pc	118.00
LMK.2620050	GC	5.0ml L-MARK Syringe FN, GT, PTFE, 23g, pt5, Headspace	pc	125.00
LMK.2620051	GC	5.0ml L-MARK Syringe FN, GT, PTFE, 26g, pt5, Headspace	pc	125.00

To order please call (800) 229-8814, (919) 929-8814 or fax (919) 929-8956

For technical support please call (877) 532-7832

www.leaptec.com

2. LC Applications (LC PAL, HTS & HTC PAL)

Part number	Market	Description	Unit	List Price
Small Volume Syringes				
PTFE(Teflon)				
LMK.2606416	LC	10ul L-MARK Syringe FN, GT, PTFE 22gs, pt 3	pc	58.00
LMK.2606416-2P	LC	10ul L-MARK Syringe, FN, GT, PTFE, 22gs	pack	109.00
LMK.2606416-6P	LC	10ul L-MARK Syringe FN, GT, PTFE 22gs, pt 3	pack	316.00
LMK.2606416C	LC	10ul L-MARK GOLD Syringe Coated FN, GT, PTFE 22gs, pt 3	pc	87.00

LMK.2620517	LC	25ul L-MARK Syringe FN, GT, PTFE, 22 gs , pt 3	pc	91.00
LMK.2620517-2P	LC	25ul L-MARK Syringe FN, GT, PTFE, 22 gs , pt 3 2/pk	pack	170.00
LMK.2620517-6P	LC	25ul L-MARK Syringe FN, GT, PTFE, 22 gs , pt 3, 6/pk	pack	496.00
100.SYRXG2522S3	LC	25ul CTC X-Type Syringe FN, GT, PTFE, 22gs, pt3, spec needle coatin	pc	205.00
LMK.2620617	LC	50ul L-MARK Syringe FN, GT, PTFE, 22 gs , pt 3	pc	91.00
LMK.2620617-2P	LC	50ul L-MARK Syringe FN, GT, PTFE, 22 gs , pt 3 2/pk	pack	170.00
LMK.2620617-6P	LC	50ul L-MARK Syringe FN, GT, PTFE, 22 gs , pt 3 6/pk	pack	496.00
LMK.2620617C	LC	50ul L-MARK GOLD Syringe Coated FN, GT, PTFE, 22 gs , pt 3	pc	138.00
100.SYRXG5022S3	LC	50ul CTC X-Type Syringe FN, GT, PTFE, 22 gs , pt3, spec needle coatir	pc	205.00
LMK.2620615	LC	50ul L-MARK Syringe FN, GT, PTFE, 22 g , pt 3	pc	91.00
LMK.2620615-2P	LC	50ul L-MARK Syringe FN, GT, PTFE, 22 g , pt 3 2/pk	pack	170.00
LMK.2620615-6P	LC	50ul L-MARK Syringe FN, GT, PTFE, 22 g , pt 3 6/pk	pack	496.00
LMK.2620615C	LC	50ul L-MARK GOLD Syringe Coated FN, GT, PTFE, 22 g , pt 3	pc	138.00
LMK.2620719	LC	100ul L-MARK Syringe FN, GT, PTFE, 22 gs , pt 3	pc	104.00
LMK.2620719-2P	LC	100ul L-MARK Syringe FN, GT, PTFE, 22 gs , pt 3 2/pk	pack	197.00
LMK.2620719-6P	LC	100ul L-MARK Syringe FN, GT, PTFE, 22 gs , pt 3 6/pk	pack	565.00
LMK.2620719C	LC	100ul L-MARK GOLD Syringe Coated FN, GT, PTFE, 22 gs , pt 3	pc	148.00
100.SYRXG10022S	LC	100ul CTC X-Type Syringe FN, GT, PTFE, 22 gs , pt3, spec needle coat	pc	205.00
100.SYRXG100223	LC	100ul CTC X-Type Syringe FN, GT, PTFE, 22 g , pt3, spec needle coatir	pc	205.00
LMK.2620721	LC	100ul L-MARK Syringe FN, GT, PTFE, 22 g , pt 3	pc	95.00
LMK.2620721-2P	LC	100ul L-MARK Syringe FN, GT, PTFE, 22 g , pt 3 2/pk	pack	179.00
LMK.2620721-6P	LC	100ul L-MARK Syringe FN, GT, PTFE, 22 g , pt 3 6/pk	pack	515.00
LMK.2620721C	LC	100ul L-MARK GOLD Syringe Coated FN, GT, PTFE, 22 g , pt 3	pc	138.00

To order please call (800) 229-8814, (919) 929-8814 or fax (919) 929-8956

For technical support please call (877) 532-7832

www.leaptec.com

2. LC Applications (LC PAL, HTS & HTC PAL) Continued

Part number	Market	Description	Unit	List Price
		PE(Polyethylene)		
LMK.2620518	LC	25ul L-MARK Syringe FN, GT, PE, 22gs, pt 3	pc	89.00
LMK.2620518-6P	LC	25ul L-MARK Syringe FN, GT, PE, 22gs, pt 3, 6/pk	pack	486.00
LMK.2620616	LC	50ul L-MARK Syringe FN, GT, PE, 22gs, pt 3	pc	93.00
LMK.2620616-2P	LC	50ul L-MARK Syringe FN, GT, PE, 22gs, pt 3 2/pk	pack	177.00
LMK.2620616-6P	LC	50ul L-MARK Syringe FN, GT, PE, 22gs, pt 3 6/pk	pack	508.00
LMK.2620616C	LC	50ul L-MARK GOLD Syringe Coated FN, GT, PE, 22gs, pt 3	pc	148.00
LMK.2620614	LC	50ul L-MARK Syringe FN, GT, PE, 22g, pt 3	pc	93.00
LMK.2620614-2P	LC	50ul L-MARK Syringe FN, GT, PE, 22g, pt 3 2/pk	pack	177.00
LMK.2620614-6P	LC	50ul L-MARK Syringe FN, GT, PE, 22g, pt 3 6/pk	pack	508.00
LMK.2620614C	LC	50ul L-MARK GOLD Syringe Coated FN, GT, PE, 22g, pt 3	pc	148.00
LMK.2620718	LC	100ul L-MARK Syringe FN, GT, PE, 22gs, pt 3	pc	104.00
LMK.2620718-2P	LC	100ul L-MARK Syringe FN, GT, PE, 22gs, pt 3 2/pk	pack	197.00
LMK.2620718-6P	LC	100ul L-MARK Syringe FN, GT, PE, 22gs, pt 3 6/pk	pack	565.00
LMK.2620718C	LC	100ul L-MARK GOLD Syringe Coated FN, GT, PE, 22gs, pt 3	pc	148.00

LMK.2620720	LC	100ul L-MARK Syringe FN, GT, PE, 22g, pt 3	pc	104.00
LMK.2620720-2P	LC	100ul L-MARK Syringe FN, GT, PE, 22g, pt 3 2/pk	pack	197.00
LMK.2620720-6P	LC	100ul L-MARK Syringe FN, GT, PE, 22g, pt 3 6/pk	pack	565.00
LMK.2620720C	LC	100ul L-MARK GOLD Syringe Coated FN, GT, PE, 22g, pt 3	pc	148.00

Large Volume Syringes

LMK.2620821	LC	250ul L-MARK Syringe FN, GT, PTFE, 22g, pt 3, barrel 6.6	pc	105.00
LMK.2620821-2P	LC	250ul L-MARK Syringe FN, GT, PTFE, 22g, pt 3, barrel 6.6 2/pk	pack	199.00
LMK.2620821-6P	LC	250ul L-MARK Syringe FN, GT, PTFE, 22g, pt 3, barrel 6.6 6/pk	pack	565.00
LMK.2620821C	LC	250ul L-MARK GOLD Syringe Coated FN, GT, PTFE, 22g, pt 3, barrel (pc	139.00
LMK.2620820	LC	250ul L-MARK Syringe FN, GT, PE, 22g, pt 3, barrel 6.6	pc	112.00
LMK.2620820-2P	LC	250ul L-MARK Syringe FN, GT, PE, 22g, pt 3, barrel 6.6 2/pk	pack	213.00
LMK.2620820-6P	LC	250ul L-MARK Syringe FN, GT, PE, 22g, pt 3, barrel 6.6 6/pk	pack	610.00
LMK.2620820C	LC	250ul L-MARK GOLD Syringe Coated FN, GT, PE, 22g, pt 3, barrel 6.6	pc	148.00
LMK.2620805	LC	250ul L-MARK Syringe FN, GT, PTFE, 22g, pt 3, barrel 8	pc	77.00
LMK.2620805-6P	LC	250ul L-MARK Syringe FN, GT, PTFE, 22, pt 3, barrel 8 6/pk	pack	425.00
LMK.2620908	LC	500ul L-MARK Syringe FN, GT, PTFE, 22g, pt 3	pc	77.00
LMK.2620908-6P	LC	500ul L-MARK Syringe FN, GT, PTFE, 22g, pt 3, 6/pk	pack	425.00
LMK.2620011	LC	1.0mL L-MARK Syringe FN, GT, 22g, pt3	pc	77.00
LMK.2620012	LC	1.0ml L-MARK Syringe FN, GT, PTFE, 23g, pt5	pc	122.00
LMK.2620054	LC	5.0ml L-MARK Syringe FN, GT, PTFE, 19g, pt5	pc	146.00

To order please call (800) 229-8814, (919) 929-8814 or fax (919) 929-8956

For technical support please call (877) 532-7832

3. Specialty Syringes

Part number	Market	Description	Unit	List Price
LMK.2620025	IFC	2.5ml L-MARK Syringe RN, PTFE, 19g, pt3	pc	146.00
LMK.2620053	IFC	5.0ml L-MARK Syringe RN, PTFE, 19g, pt3	pc	135.00
LMK.2822412	IFC	Removable needle for 2.5ml & 5.0ml Syringes, 19g, pt 3	pc	28.00
LMK.2620510	GC LINEX	25µl L-MARK Syringe ø6,5 FN 0,47(G26s)c70	pc	
	DMI			93.00
	Collect		pc	
LMK.2620827	PAL	250ul L-MARK Syringe f/ Collect PAL with 10-32 fittings at inlet and outl		275.00
Replacement Plungers				
LMK.2636402	LC & GC	10ul L-MARK replacement plunger, PTFE tip	pc	23.00
LMK.2636605	LC & GC	50ul L-MARK replacement plunger, PTFE tip	pc	45.00
LMK.2636605-2P	LC & GC	50ul L-MARK replacement plunger, PTFE tip 2/pk	pack	84.00
LMK.2636605-6P	LC & GC	50ul L-MARK replacement plunger, PTFE tip 6/pk	pack	245.00
LMK.2636708	LC & GC	100ul L-MARK replacement plunger, PE tip	pc	53.00
LMK.2636708-2P	LC & GC	100ul L-MARK replacement plunger, PE tip 2/pk	pack	102.00
LMK.2636708-6P	LC & GC	100ul L-MARK replacement plunger, PE tip 6/pk	pack	290.00
LMK.2636606	LC	50ul L-MARK replacement plunger, PE tip	pc	52.00
LMK.2636606-2P	LC	50ul L-MARK replacement plunger, PE tip 2/pk	pack	96.00
LMK.2636606-6P	LC	50ul L-MARK replacement plunger, PE tip 6/pk	pack	284.00
LMK.2636707	LC	100ul L-MARK replacement plunger, PTFE tip	pc	59.00
LMK.2636707-2P	LC	100ul L-MARK replacement plunger, PTFE tip 2/pk	pack	112.00
LMK.2636707-6P	LC	100ul L-MARK replacement plunger, PTFE tip 6/pk	pack	323.00
PAL.PLGX25	LC	25uL CTC X-Type Replacement Plunger for X-Type syringe, 3/pk.	pack	225.00
PAL.PLGX50	LC	50uL CTC X-Type Replacement Plunger for X-Type syringe, 3/pk.	pack	225.00

PAL.PLGX100 LC 100uL CTC X-Type Replacement Plunger for X-Type syringe, 3/pk. pack 225.00

LEGEND

NOTE: if not noted, needle length is 51mm

AS Cone shape needle tip

pt3 Straight cut for LC

AS Cone tip for GC

pt5 Pointed needle tip with side hole

FN Fixed needle

RN Removable needle

GT Gastight

PTFE Teflon® (Polytetrafluoroethylene)

PE Black Polyethylene

Barrel Diameter Barrel diameter- small syringes (5-100ul) are in general 6.5mm; larger syringes (250ul and 500ul 7.7mm)

Needle normal ID-g

Gauge small ID-gs or s

To order please call (800) 229-8814, (919) 929-8814 or fax (919) 929-8956

For technical support please call (877) 532-7832

www.leaptec.com

® Registered Trademark of LEAP Technologies, Inc. All rights reserved.