

# Nano PAL

(not expandable)



- LC PAL w/ single valve drive
- Single micro-well cooled tray holder
- Run integrated with Analyst™, Eksigent NanoLC or Flux Rheos 2000

# Nano PAL Specs

- Single position cooled tray holder for 1 deep well plate, shallow 96 well plate or 98 2ml vials



- Fast wash station for syringe washes with 2 separate solvents



- LEAP Shell for integration with Analyst™ and Flux Rheos 2000 with splitter and mass flow monitor



- PEAK micro valve for low carry over injections to Nano Spray MS

